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

## 1.1 Key

### WEINMANNstrap


- 1 Forehead cushion
- 2 Mask seal
- 3 Head strap
- 4 Rotating sleeve
- 5 Support
- 6 Seal for ports
- 7 Mask body
- 8 Elbow
- 9 Exhalation connector
- 10 Exhalation slot
- 11 Hose connector

### SOMNOstrap

### HEADstrap

## 1.2 Purpose of Use

The nasal mask with exhalation connector is an accessory that serves as a sealing element between patient and therapy device. It is used in the treatment of sleep apnea and for non-invasive ventilation of patients with respiratory insufficiency.

## 1.3 Functional Description

The nasal mask and the head strap are attached over your nose so that the flow of air from the therapy device is routed through the nasal mask to the nose. The form and flexible material of the mask seal creates an airtight bond to the face. The elbow connection allows the greatest possible freedom in positioning the hose.

The inhaled air flows through the mask/hose system. The exhaled air escapes through a slot between the elbow and the rotating sleeve. This allows you to breathe quietly during therapy. The exhalation connector avoids the risk of re-inhalation.

## 2 For Your Safety

### Safety regulations

For your and your patients' safety and in accordance with the requirements of the European Union Guidelines 93/42/EEC, please observe the following points:

- Please read these instructions for use very carefully. They are a part of the mask and should be available at all times.
- If third-party items are used, functional failures may occur and fitness for use may be restricted. Biocompatibility requirements may also not be met. Please note that in such cases, any claim under warranty and liability will be voided if neither the accessories nor genuine replacement parts recommended in the instructions for use are used.
- The mask is designed for therapy pressure between 4 and 30 hPa. The nasal mask must not be used at higher or lower therapy pressures.
- Note that patients with limited spontaneous respiration must be monitored at all times and that the appropriate negative pressure/leak alarm on the ventilation device must be activated.
- Patients with obstructive and restrictive lung diseases should be examined on an individual basis during use of the *SOMNOplus*, for example with regard to re-inhalation of CO<sub>2</sub>.
- The *SOMNOplus* has an integrated exhalation system. Do not connect any other exhalation system or exhalation valve to the nasal mask.
- Do not block any of the openings in the nasal mask.
- Only put the nasal mask over your face for a prolonged period when the therapy device is switched on. The flow of air from the therapy device serves among other things to remove the exhaled air from mask and hose. If the therapy device is not switched on, there is a risk of re-inhalation when you breathe through your nose.
- Do not use a breathing hose that is electrically conductive (exception: breathing tubes which can be heated electrically).

- Clean the mask as described under Section “**4 Hygienic preparation**”. Inappropriate cleaning procedures, unsuitable cleaning agents and infrequent cleaning can result in an infection or bacterial contamination and can impede the functioning of the mask.

## Side effects

The following side effects may occur during short-term or long term use:


nasal congestion, dry nose, dry mouth in the morning, sinus pressure, irritation of the conjunctiva as well as irritating noises during exhalation. Wearing the mask can also lead to tenderness of the face or reddening of the skin.

# 3 Operation

## Assemble and adjust forehead support


The forehead support is made up of a support and a forehead cushion. The support makes it possible to adjust the mask to your face while the forehead cushion allows for comfortable distribution of the mask pressure across the forehead. By adjusting the support, pressure on the bridge of the nose can be minimized (Fig. **A**).

**A**

Select the support size that best suits your face (Fig. **B**). Attach the forehead cushion to the support and place the support in the slit on the mask body. In most cases, the size 3 support will be appropriate. The forehead cushion is already set up on this support when you receive the mask.

**B**


Place the mask on your face and check to see if you have chosen the correct support.

If the mask isn't tight or if it causes uncomfortable tenderness, shift the support within the slot of the mask or choose a different support size.

### Dismantle forehead support

Loosen the forehead support by pressing on the clamp on the mask body (Fig. C). Pull forehead cushion off support.

C


### Head strap adjustment

When you first put on the nasal mask with the head straps, it is a good idea to ask another person for help.


- Thread the four Velcro straps of the head strap into the slots on the mask body. You can use either the lateral slots or the top hoop to attach the Velcro forehead straps.
- Pass the head strap over your head and hold the nasal mask over your nose.
- Adjust the head strap with the aid of the Velcro straps so that the mask seals on your face but is not too tight (Fig. D1).

D1

- To adjust the Velcro strap on top of the head: release the Velcro strap, pull gently on the ends of the head strap with both hands and push the Velcro strap back on (Fig. D2).


D2

## Connection to therapy device

Push the hose of your therapy device onto the rotating sleeve (Fig. E). In doing so, ensure that there is a secure, airtight connection to the hose. The connection should not come loose as a result of therapy pressure.

**E**

If you wish to interrupt therapy briefly, you can pull the rotating sleeve out of the elbow. It is not necessary to remove the mask or pull the breathing hose off the rotating sleeve with this nasal mask.

## Combining with other devices

This nasal mask can be combined with any device to perform the therapies described under “**1.2 Purpose of Use**”, as long as the following conditions are met:

- when combining the nasal mask with other devices, be aware that the therapy pressure in the nasal mask depends on the combination of devices. Ask your doctor or dealer whether the correct therapy pressure can be achieved in the mask using the desired combination of devices. Have the settings of your therapy device checked by your doctor or dealer.
- take account of the specifications in the instructions for use for the devices you intend combining with the nasal mask.


## Use with cold-air or warm-air humidifier

The nasal mask system can be used with a cold-air or warm-air humidifier. The degree of humidity you choose should not be so high that water condenses in the air hose. Please be sure to read the instructions for use for your humidifier and your therapy device.

## Pressure measurement

SOMNOplus is supplied with two ports. They can be used for external pressure measurement.

If only one port is required, turn the connection upwards (Fig. F). When the functions with the ports are no longer required, they should be made airtight with the enclosed seal. Otherwise air can escape through the ports and therapy pressure will not be maintained.


**F**

### After use

- Undo one Velcro strap on the head strap and remove the mask.
- Clean the mask as described in Section "4 Hygienic preparation".

### Dismantling the mask

1. Remove the head strap from the nasal mask.
2. Loosen the forehead support from the mask body by pressing on the clamp. Pull the forehead cushion off the support.
3. Pull the mask seal off the mask body and remove the seal for the ports.
4. Pull the rotating sleeve out of the elbow and remove the elbow from the mask body by pressing gently inwards (Fig. **G**).

**G**

### Assembling the mask

1. Push the elbow into the mask body from inside. Push the rotating sleeve into the elbow. Make sure both parts engage properly.
2. Push the mask seal onto the mask body and the port seal onto the pressure measuring port.
3. Attach the forehead support to the mask body as described under "Assemble and adjust forehead support".
4. Attach the head strap to the nasal mask as described under "Head strap adjustment".

## 4 Hygienic preparation

### 4.1 Schedule

#### Daily cleaning

Please clean the nasal mask after every usage as described in item **"4.2 Procedure for cleaning, disinfecting or sterilizing"** under "Cleaning by hand". For cleaning purposes we recommend taking the mask apart as described above.

#### Weekly cleaning, disinfecting or sterilizing

The nasal mask should be cleaned thoroughly at least once per week.

For cleaning purposes, take the mask apart as described above.

### 4.2 Procedure for cleaning, disinfecting or sterilizing

The following procedures are allowed:

Procedures					
Parts	Cleaning by hand	Dishwasher	Washing machine	Disinfecting*	Sterilizing (as an alternative to disinfecting))
Mask seal	X	X		X	X
Forehead cushion	X	X		X	X
Seal for ports	X	X		X	X
Mask body	X	X		X	
Support	X	X		X	
Elbow	X	X		X	

Procedures					
Parts	Cleaning by hand	Dishwasher	Washing machine	Disinfecting*	Sterilizing (as an alternative to disinfecting))
Rotating sleeve	X	X		X	
Head strap:					
- WEINMANNstrap	X				
- HEADstrap	X				
- SOMNOstrap	X		X		

Cleaning by hand	Clean in warm water, free of residues with a mild detergent as for dishwashing. Carefully rinse the parts with clear water.
Dishwasher	Wash at max. 65°C. Carefully rinse the parts with clear water.
Washing machine	Wash at 30°C.
Disinfection at home (no change in patients)*	Clean all parts as described under "Cleaning by hand" and then immerse them in a diluted solution of disinfectant. After disinfecting rinse all parts thoroughly with clear water.
Disinfecting at clinics and with change in patients*	Clean in a diluted solution using mechanical means, e.g. with a brush. After disinfecting, rinse the parts thoroughly with distilled water and dry them completely before you reassemble the mask.
Sterilizing (as an alternative to disinfecting)	Steam-sterilize in devices complying with EN 285. Temperature 134 °C, Minimum time: 3 minutes

Dry all parts thoroughly after cleaning. Do not dry parts in direct sunlight. Wash the headgear before using for the first time, as the colour might run. The headgear may only be washed by hand! Do not iron the headgear, otherwise the hook-and-loop closures will no longer work! Do not dry the headgear in a dryer. Check components for residual dirt or damage. Repeat cleaning if necessary. If parts are damaged, replace them. Assemble the nasal mask as described above.

**\* We recommend the disinfectant gigasept® FF (new). Make sure that you follow the manufacturer's directions. Tolerance of this disinfectant is proven. We recommend that you wear suitable gloves when disinfecting the equipment (e.g. household or disposable gloves).**

Notes on the hygiene treatment and on the hygiene treatment for a change of patient can be found in a brochure from our information portal at [www.weinmann.de](http://www.weinmann.de). If you do not have Internet access, contact Weinmann.

### 4.3 Instructions for Disinfecting and Sterilizing

For patients with a weakened immune system or a special medical history daily disinfection or sterilisation may be necessary as for a change in patients. In this case, please consult your doctor.

## 4.4 Disinfecting and Sterilizing for new patient

If the nasal mask is to be used for another patient, it must be subjected to hygiene measures first. Dismantle the mask and disinfect or sterilize all parts as described above or use a new mask. Use a new head strap.

## 5 Functional checks

Please carry out a visual check once per week. If the mask seal is no longer airtight due to a tear, please change it.

Please check the functional capacity and airtightness of the ball-and-socket joint. If it no longer moves easily or is no longer airtight, please change the elbow connection and/or the body of the mask.

Parts that can no longer be used can be disposed of normal household waste.

## 6 Troubleshooting

Fault	Cause	Remedy
In the morning you have pain from pressure on your face.	The nasal mask is too tight.	Enlarge the headgear or the head straps a bit. Change the position of the forehead support on the mask body. If necessary, select a different support size for the forehead cushion.

Your eyes feel a draught.	The mask is too loose	Make the headgear or the straps a little tighter. Select a different support size for the forehead cushion.
	The mask doesn't fit.	Please see if you need a different size mask. Your specialist dealer can advise you.
Therapy pressure is not reached.	The mask is not airtight.	Correct the fit, adjust the head straps. Change the position of the forehead support on the mask body. If necessary, select a different support size for the forehead cushion.
	Tears in the shaped mask.	Change the mask seal.
	Supply system is not airtight.	Check the plug connection and the fit of the hose muff.
	Mask blows air over the open ports.	Place the seal on the ports.

## 7 Supply Schedule

	with <b>SOMNOstrap</b>	with <b>WEINMANNstrap/ HEADstrap</b>
SOMNOplus, size S consisting of: - Mask seal size S - Body of mask size S - Seal for ports - Forehead support sizes 1-4 (size 3 fitted) - Exhalation connector, consisting of elbow and rotating sleeve - EASYclip - WEINMANNstrap  - SOMNOstrap blue  - HEADstrap - Instructions for use	WM 25710	WM 25710/1
	WM 26210	
	WM 25112	
	WM 25165	
	WM 25180	
	WM 25705	
	WM 15699	
		WM 26367
	WM 23695	
		WM 26351
SOMNOplus, size M consisting of: - Mask seal size M - Body of mask size M	WM 25720	WM 25720/1
	WM 26220	
	WM 25122	
SOMNOplus, size L consisting of: - Mask seal size L - Body of mask size L	WM 25730	WM 25730/1
	WM 26230	
	WM 25132	

### Replacement parts

You can order accessories and spare parts separately if required. A current list of accessories and spare parts can be ordered on the Internet at [www.weinmann.de](http://www.weinmann.de) or via your dealer.

WM 16806f 09/2012

## 8 Technical Data

Device class acc. to 93/42 EEC:	II a
Dimensions (WxHxD):	approx. 80 x 140 x 90 mm
Weight:	approx. 55 g
Therapy pressure range:	4 to 30 hPa
Hose connector:	Cone 22 mm diam. (EN ISO 5356-1) for elastomeric hoses
Port pressure measurement:	Ø 4 mm
Connection head straps:	Strap width max. 28 mm
Operating temperature range:	+ 5 °C to + 40 °C
Storage temperature range:	- 20 °C to + 70 °C
Flow resistance at 50 l/min	0.06 hPa
Flow resistance at 100 l/min	0.22 hPa
Quoted two-number noise emission value to ISO 4871: - sound pressure level - sound power level (uncertainty factor: 3 dB(A))	15 dB(A) 23 dB(A)
Duration of use	6 to 12 month depending on frequency of use and cleaning

**CE 0197**

Subject to technical change without notice


## Materials

Mask seal	silicone
Mask body	polycarbonate
Elbow	polycarbonate
Rotating sleeve	polycarbonate
Seal for ports	silicone
Forehead cushion	silicone
Support	polypropylene
WEINMANNstrap	Elasthane, Polyester, Polyurethane, Nylon, cotton, polyamide
SOMNOstrap	cotton, polyamide
HEADstrap	Nylon/Spandex, Polyurethane, Polyester/ Elasthane

All individual parts of the nasal mask are latex-free.

## Characteristic pressure/flow graph

The characteristic pressure/flow graph shows outlet flow as a function of therapy pressure.

## 9 Warranty

Weinmann gives the customer a limited manufacturer warranty on new genuine Weinmann products and any replacement part fitted by Weinmann in accordance with the warranty conditions applicable to the product in question and in accordance with the warranty periods from date of purchase as listed below. The warranty conditions can be downloaded from [www.weinmann.de](http://www.weinmann.de) on the Internet. We can also send you the warranty conditions on request.

In the event of a claim under warranty, contact your specialist dealer.

<b>Product</b>	<b>Warranty period</b>
Masks including accessories, rechargeable batteries, batteries (unless quoted differently in the technical documentation), sensors, tube systems	6 months
Disposable products	None

## 10 Declaration of conformity

Weinmann Geräte für Medizin GmbH + Co. KG declares herewith that the product complies fully with the respective regulations of the Medical Device Directive 93/42/EEC. The unabridged text of the Declaration of Conformity can be found on our website at: [www.weinmann.de](http://www.weinmann.de)

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